

100 bp DNA Ladder , ready-to-use

Description:

100 bp DNA Ladder ready to use is used for determining the size of double-stranded DNA from 100 to 1500 base pairs on agarose gels. The 500 bp fragment is present at increased intensity to allow easy identification. 100 bp DNA Ladder is supplied in storage buffer with additional dye for easy and direct application to the gel. 5x DNA Loading Buffer can be used to add dye to the testing samples in a final concentration of 1x.


Number of fragments: 12 (blunt-end), stored in 1x Loading Buffer

DNA fragments length (bp):

1500, 1200, 1000, 900, 800, 700, 600, 500, 400, 300, 200, 100

Usage recommendation: 5 µl per lane

Storage: -18 to -22 °C (long-term)
+2 to +4 °C (up to 6 months)



REF	306120	colour
100 bp DNA Ladder ready-to-use	500 µl	transparent
5x DNA Loading Buffer	1 ml	transparent

Concentration gradients may form in frozen product. Mix well before using the DNA ladder.

Do not heat the ladder above 37 °C!

BIORON GmbH – In den Rauhweiden 20 - 67354 Römerberg (Germany)

Phone +49 6232 298 45 0 - Fax +49 6232 298 45 29

info@bioron.net www.bioron.net